

KN95 FILTER MASK



Description

KN95 Filter Mask 10 Pack
Standards: GB-2626-2006
Certificate: CE EN149:2001+A1:2009
Filter Performance: $\geq 95\%$

Materials

Number of Layers: 5
Shell: Non-woven Fabric
Filter: Meltblown (x2)
Inner: Non-woven Fabric
Ear Loop: Spandex
Nose Clip: Aluminum

Others

Dimensions (Folded) 155x105mm (6.1"x 4.1")
Weight 5.3g (0.18 oz)

Circadiance offers this mask to help our customers meet your personal protective equipment needs. Contact us for details.

Orders will be accepted for 10 Packs (Part #101738) at orders@circadiance.com or +1 888 825 9640.

