

# DISPOSABLE MASK



## Description

Disposable Mask 50 Pack  
Standards: YY/T 0969-2013  
Certificate: CE EN149:2001+A1:2009  
BFE:  $\geq 95\%$

## Materials

Number of Layers: 3  
Shell: Non-woven Fabric  
Filter: Melt blown  
Inner: Non-woven Fabric  
Ear Loop: Spandex  
Nose Clip: Aluminum

## Specs

Dimensions: 175\*95mm (6.9" x 3.7")  
Weight: 3.2g (0.11 oz.)

**Circadiance offers this mask to help  
our customers meet your personal  
protective equipment needs.  
Contact us for details.**

**Orders will be accepted for 50 Packs (Part #101746)  
at [orders@circadiance.com](mailto:orders@circadiance.com) or +1 888 825 9640.**

